FR

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The Instructions for Use (IFU) is related to ALCIS clamp from range 2023-PIN. It is an additional part of the IFU associated to ALCIS depth electrodes from the range 2069-EPC and associated accessories (ranges 2023 and 2036).

1. Products description

Technical characteristics of Clamp:

2023-PIN Year of CE marking : 2009	Clamp	  <i>Technical form 2023-AVG-C-FT</i>	Use exclusively with ALCIS short guidance screws range 2023-VG-C .
---	-------	---	--

Assigned performances

Ability to maintain the short guidance screw in place behind the stereotactic guide while screwing.

2. Intended use

Clamp is intended to maintain the Alcis short guidance screw in place behind the stereotactic grid when screwing. The target population is drug resistant epilepsies patients.

3. Contraindications

No specific contraindications.

4. Recommendations / Warnings

- Do not "force" for accessories assembly or disassembly (assemblies must be performed without excessive stress)
- The device is intended to be used by experienced neurosurgeons.
- Check the integrity and the condition before each use.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state in which the user/the patient is established.
- Be careful to not drop the device because it shall result in the deterioration/breakage of the clamp.

Storage and maintenance

- Store in a cold dry place and away from sunshine at room temperature.
- Inspect packaging before use to detect any damage.
- Do not use in case of wear and/or damage of any kind.

Cleaning, Disinfection and sterilization

Reusable accessories

any devices identified with any of the symbols shown here

(Refer to product label).

- Reusable devices are provided unsterile
- Before use, clean and sterilize according to a controlled and validated process (alkaline detergent and sterilization in an autoclave at 134°C for 18 minutes – prion cycle).
- During the cleaning take care to biological residues absence.
- After use, clean the reusable ancillaries (alkaline detergent).

Disposal procedure

There are infectious risks related to clamp disposal, according to the WHO waste classification. To minimize the risks, the most common treatments for infectious waste are: steam treatment (autoclaving), integrated steam treatment (hybrid autoclaves or advanced steam treatment), microwave treatment, dry heat treatment, chemical treatment (internal shredding or chemical disinfectants) and incineration (combustion, pyrolysis or gasification). One of this method shall be followed to reduce the infectious risk.

5. Main complications

No specific complications or potential adverse effects are identified for the clamp.


Main complications of interventional technics in stereotactic conditions still apply.

6. Warranty limit

ALCIS products are submitted to control tests at several stages of manufacturing process to ensure compliance with specifications. ALCIS can never be held responsible for any expenses or damages directly or indirectly resulting from:

- The non-respect of instructions, recommendations and warnings of the instruction for use,
- The use of accessories or devices unsuited or defective

7. Use procedure

D_a being the distance between the proximal side of the guide (grid) and the skin, and D_b the distance between the extremity of the drill and the stopper, determine L_s the length of guidance screw to use by choosing the closer length to the bone + skin thickness ($D_b - D_a$).

Hold the short guidance screw with the clamp **2023-PIN** at the entrance hole.

Place the long screwdriver 2023-TO-C-160 through the grid or the guide and insert it into the guidance screw (its hexagonal part entirely entered into the screw). Without the clamp, fully screw in the osseous hole (screw is self-threading).

References table

Reference	Length L (mm)	Material	Packaging
2023-PIN	135	Stainless steel	Non sterile units

8. Definition of labelling and packaging symbols

REF	Manufacturer's Catalog Reference	Keep away from sunlight	LOT	Batch No
	Store in a dry place	Do not use if damaged	x50	Sterilize to re-use validated for x <i>nn</i> cycles (here 50)
	Manufacturing date	24M		Read Instruction For Use
	Manufacturer			Warning Steam sterilize before use: 134°C – 18min
	International phone number			Medical Device

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