FR

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
EN

This Instructions for Use (IFU) is related to ALCIS screwdriver from range 2023-TO. It is an additional part of the IFU associated to ALCIS depth electrodes from the range 2069-EPC and associated accessories (ranges 2023 and 2036).

1. Product description

Technical characteristics of Screwdriver:			
Range 2023-TO  Year of CE marking : 2002	Screwdriver 	  REF 2023-TO-C-III <i>III being the length</i> <i>Technical form 2023-AVG-C-FT</i>	Body OD adapted to ID 2.5mm guides To be used exclusively with ALCIS short guidance screws range 2023-VG-C . Available in length 20 or 160mm
		REF 2023-TO-III <i>III being the length</i> <i>Technical form 2023-APR-FT</i>	To be used exclusively with ALCIS standard guidance screws range 2023-VG . Available in length 20 or 160mm
		  REF 2023-TO-BC <i>Technical form 2023-AVG-C-FT</i>	To be used exclusively with ALCIS short guidance screws range 2023-VG-C . Available in length 30mm
Range 2023-TO-BC  Year of CE marking : 2002	Screwdriver for short guidance screw's cap 		

Performances:

- To be assembled to Alcis guidance screw in order to screw or unscrew them on the patient skull and, if applicable, through a guide with compatible diameter.
- To be extracted from the guidance screws without any particular difficulty.
- For 2023-TO-BC: to screw/unscrew the cap of the short guidance screw without any particular difficulty.

2. Intended use

The screwdriver is intended to screw and unscrew previously unlocked Alcis guidance screw. The screwdriver for short guidance screw's cap is intended to screw/unscrew the cap in case of difficulty to remove it.

As the screwdrivers are accessories of Alcis Depth Coagulation Electrodes and guidance screws, the target population is drug-resistant epilepsies patients.

The intended users are experienced neurosurgeons.

3. Contraindications

No specific contraindications.

4. Recommendations / Warnings

Warning

- Do not "force" for accessories assembly or disassembly (assemblies must be performed without excessive stress)
- Check accessories integrity before each use.

Storage and maintenance


- Store in a dry environment at room temperature.
- Inspect packaging before use to detect any damage.
- Do not use in case of wear and/or damage of any kind.

- Do not use if package is opened or damaged.

Cleaning, disinfection and sterilization

Reusable accessories:

any devices identified with any of the symbols shown here

(Refer to product label).

- Reusable devices are provided unsterile
- Before use, clean and sterilize according to a controlled and validated process (alkaline detergent and sterilization in an autoclave at 134°C for 18 minutes – prion cycle).
- During the cleaning, take care to biological residues absence.
- After each use, clean the reusable ancillaries (alkaline detergent).

any devices identified with symbol shown here (refer to product label)

Disposal procedure :

There are infectious risks related to screwdriver disposal, according to WHO waste classification. To minimize the risks, the most common treatments for infectious waste are: steam treatment (autoclaving), integrated steam treatment (hybrid autoclaves or advanced steam treatment), microwave treatment, dry heat treatment, chemical treatment (internal shredding or chemical disinfectants) and incineration (combustion, pyrolysis or gasification). One of this method shall be followed to reduce the infectious risk.

5. Main complications

Main complications of interventional technics in stereotactic conditions still apply.

Note: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state in which the user/the patient is established.


6. Warranty limit

ALCIS products are submitted to control tests at several stages of manufacturing process to ensure compliance with specifications. ALCIS can never be held responsible for any expenses or damages directly or indirectly resulting from:

- The non-respect of instructions, recommendations and warnings of the instruction for use,
- The purchase, the use, the implantation, the removal or replacement of screwdriver,
- The use of accessories or devices unsuited or defective

7. Use procedure

During the positioning of guidance screws

- D_a : distance between the proximal side of the guide (grid) and the skin.

- D_b / distance between the extremity of the drill and the stopper.


- Determine L_s the length of guidance screw to use by adding **10mm** to the bone + skin thickness ($D_b - D_a + 10mm$).

- Insert the guidance screw by square end into a long screwdriver.


- Insert this assembly through the grid or guide in the osseous hole *screw is self-threading*.

- **Screw the screw.**

Using short guidance screw (2023-VG-C)


Determine L_s the length of guidance screw to use by choosing the closer length to the bone + skin thickness ($D_b - D_a$).

Hold the short guidance screw with the clamp 2023-PIN at the entrance hole.

Place the long screwdriver through the grid or the guide and insert it into the guidance screw (its hexagonal part entirely enter into the screw). Without the clamp, fully screw in the osseous hole (screw is self-threading).

- Secure screwing depth with a suitable depth-stopper placed at a distance equal to D_b from the distal end of guidance screw.

- Reiterate the above operations for each screw.

We advise to control the screws position (matching with planned trajectories) with a suitable technique.

During electrodes implantation

- Insert the electrode through the guidance screw up to the seal, then screw tight the cap.

Pay attention to avoid to move the sliding seal during this step.

- Check with a slight traction that the electrode is well fixed.

Using short guidance screw (2023-VG-C)

Adjust the seal location to $L_i + 2,5$ mm, L_i being the implanted length


Mount the cap (provided in the package) on the electrode.

Then insert the electrode, fully screw the cap maintaining the sliding seal close to the cap (if necessary use the screwdriver 2023-TO-BC adapted for the cap).

Pay attention to avoid to move the sliding seal during this step.

Check with a slight traction that the electrode is well fixed.

Maintain the electrode in a groove of the cap and put the ring to keep the electrode blocked in position.

During electrodes removal

Warning: To do without delay at the end of the thermo-coagulation procedure

- Check the absence of abnormalities with a suitable technique.
- Unscrew the cap, then pull the electrode to remove it.

Using short guidance screw (2023-VG-C)

Remove the retaining ring and unscrew the cap using a suitable screwdriver (**2023-TO-BC**)


Pull the electrode to remove it

- Unlock the screw with a wrench (2023-CL-C), then totally **unscrew** with a short screwdriver.
- Close the cutaneous hole.
- Reiterate the same operations for each electrode.
- Check the cerebral structures with a suitable technique.

References table

Reference	Material	Packaging
2023-TO-C-020	Stainless steel (rod) POM (handle)	Non sterile units
2023-TO-C-160		
2023-TO-020		
2023-TO-160		
2023-TO-BC	Stainless steel	Non sterile units

8. Definition of labelling and packaging symbols

Table 1 – Symbols definition				
 REF	Manufacturer's Catalog Reference	 Keep away from sunlight	 LOT	Batch No
	Store in a dry place	 Do not use if damaged		Sterilize to re-use validated for x <i>nn</i> cycles (here 50)
	Manufacturing date	 24M		Read Instruction For Use
	Manufacturer	 Warning		Steam sterilize before use: 134°C – 18min
	International phone number	 MD		Medical device

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