FR

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DEPTH COAGULATION ELECTRODES

Range 2069 – EPC

The recommendations of the Instructions For Use (IFU) come from the SEEG state-of-the-art, as well as experiments with Alcis depth coagulation electrode - range 2069-EPC. They apply to all products listed in the Table 2 at the end of this document.

The use of ALCIS depth electrodes needs to take note of the following elements:

	Depth electrodes IFU	Ref. 2069-EPC-001
	Foramen oval introducer IFU	Ref. 2069-EPC-005
Implantation, recording & removal	Ruler IFU	Ref. 2023-TG-013
	Wrench IFU	Ref. 2023-CL-IFU
	Screwdriver IFU	Ref. 2023-TO-IFU
	Clamp IFU	Ref. 2023-PIN-IFU
Coagulation	Instructions for coagulation IFU	Ref. 2069-EPC-IFC
	Connecting box IFU	Ref. 2069-BC3-IFU

1. Products description

Intended use

The ALCIS depth electrodes (range 2069-EPC), are intracerebral multicontacts electrodes intended for pre-surgical diagnosis of drug-resistant epilepsies. They are designed to:

- **stereo-electroencephalography recording (SEEG), monitoring, and short stimulation**, in association with recording, monitoring and stimulation systems,
- **coagulation « in situ »**, in association with RF generator, as a therapeutic alternative to conventional epilepsy surgery.

***Check with your local reseller if coagulation application is legally registered in your country.
For more information refer to the instruction for coagulation (ref. 2069-EPC-IFC).***

Performances

Different types of contacts are available: standard macro, short macro and micro contacts. The number, type and distribution of the contacts vary according to the model to allow a choice adapted to the brain zone to explore. The implantation period vary between 1 to 2 weeks. Short terms use limited to less than 30 days.

The depth electrodes are designed to be implanted in surgical context with standard stereotactic means or by foramen oval way, to collect and transmit in association with an extension cable, the electrical signal of brain activity to a recorder system or from a stimulator to the brain with a minimal alteration. They can be spotted by classical XRay imaging means. They can be fixed on skull with guidance screw for an effective fixation and tightening. They allow, in association with compatible accessories, to thermo-coagulate *in situ* between two macro contacts the surrounding tissues.

• Stopper Ref 2036-BU-XX-RU

- Ability to be assembled with a compatible diameter to an adjustable distance (\approx desired cut depth) in order to ensure a stopping effect.
- Ability to ensure a minimum resistance of 20N when screwed at abutment.
- Ability to be steam sterilized (134°C – 18 minutes)

• Extension Cable Ref 2069-CE


- Ability to connect and convey electrical signals with minimal alteration, between the recorder / stimulator and the depth coagulation electrode
- Ability to be used and re-used in the best safe conditions for patients and users

Contents of the packaging

Any 2069-EPC products include:

- this Instruction For Use (is enclosed in each packaging);
- traceability labels (enclosed in the box label)
- an ALCIS depth electrode (provided with a repositionable seal and a cap for depth adjusting, tightening, and fixing on guidance screw).

Kits packaging includes a guidance screw or a foramen oval introducer as illustrated below.

2. Indications

The use of ALCIS depth electrodes is indicated for:

- **The recording and/or the short electrical stimulation**, for the pre-surgical diagnosis of drug-resistant epilepsies.
- **The bipolar thermo-coagulation *in situ*** between two contiguous macro contacts, (see instruction for coagulation), for the therapeutic treatment of epilepsy, especially for patients for whom conventional surgery is excluded. (For the time being, thermocoagulation with this type of electrodes was made in cases of cortical dysplasia, sclerosis of the hippocampus, periventricular heterotopia, cryptogenic epilepsy and postmeningitic epilepsy).

Check with your local reseller if coagulation application is legally registered in your country.

3. Contraindications

- Any other use of depth electrodes 2069-EPC than those described in this document is contraindicated.
- The contraindications related to allergies to electrodes materials are applicable.
 - See § 4 Materials in contact with the patient.
- The contraindications related to neurosurgical interventions techniques in stereotactic conditions are applicable.
 - The neurosurgeon must verify the adequacy of pathological and anatomical condition of the patient with the recording, stimulation and coagulation procedures via ALCIS depth electrodes.
- The use of other electrotherapy technics than RF thermo-coagulation is formally contraindicated on a patient carrying neurological electrodes.
 - The use of other techniques using the electromagnetic waves on a patient carrying neurological electrodes must be submitted to the agreement of the doctor responsible for the implantation and/or of the patient follow-up (AFSSAPS* note, 2001/06/18).

*French health authority in 2001.

4. Recommendations / Warnings

The implantation, the stimulation and the coagulation with ALCIS electrodes must be carried out by an experienced neurosurgeon and according to a controlled stereotactic procedure. The choice of the model and the number of electrodes to be established, as well as the type of coagulations to create are responsibility for the neurosurgeon. The patient must be placed under a continuous surveillance and in an adapted environment to avoid every risk of incidents.

- Handle Depth Coagulation Electrodes and any associated accessories with care.
- No modification of any device or equipment covered by this IFU is allowed.
- Check that the equipment to which the coagulation electrode 2023-EC is connected are operational and meet the international standards of safety. Notably: HF Generator used must comply with IEC/EN 60601-1 & IEC/EN 60601-2-2 and associated applicable standards).
- Check that the equipment to which the depth electrode 2069-EPC are connected are operational and meet the international standards of safety. Notably: SEEG Recorder used must comply with IEC/EN 60601-1 & IEC/EN 60601-2-26 (and associated applicable standards).

Implantation, Recording & Stimulation

- Plan electrodes trajectories in vessel-free corridors.
- Ensure that an available product is always ready for replacement in case of dysfunction of the first used.
- Do not twist the cables on themselves.
- Check the integrity and the good state of electrodes and all accessories before use (including fixing accessories).
- Do not use a product in case of wear and / or damage of any kind.
- Before the first implantation check the compatibility of recording and short stimulation systems with the ALCIS depth electrodes and extension cables. It is the responsibility of the user to ensure himself before use of the mutual compatibility of the products.
- It is the responsibility of the user to use equipment complying with the local regulation.
- To avoid any risk of deterioration of electrode performances, **always respect the minimal length of implantation indicated on label and on product technical form.**
 - skin or dura coagulation with 2023-EC:
 - ➔ Maximal tension allowed 2kV
 - ➔ 30W advised in coagulation mode
 - It is recommended to use antistatic gloves.
 - It is not recommended to use needle monitoring electrodes.
 - Implanted, the depth electrode establishes a direct conductive connection with patient.
 - ➔ Avoid any contact with any powered device or any device that could be powered.
 - ➔ Keep as far as possible any powered device or devices that could be powered.

After linking to any electro medical system, the distal part of electrode becomes an applied part type BF.

Removal

- Explant gently the electrodes.
- Guidance screw: avoid starting unscrewing with screwdriver. Use a suitable wrench to start.

Incompatibility with other techniques

- **See §3 Contraindications.**
- During the MRI imaging, to reduce the risks of induced current in the device:
 - ➔ Avoid forming loops with the not-implanted part of the electrode,
 - ➔ Disconnect the extension cable and the electrode.
- To date, the absence of burns risks during MRI imaging haven't been clearly establish. So the decision to carry out MRI or MEG imaging to a patient with neurological electrode is under neurosurgeon responsibility. We advise to use Scanner or Xray imaging.
- Constitutive materials of electrodes and of fixing system are adapted to an MRI or MEG examination. Nevertheless they may create artefacts of imagery which have to be taken into account in the interpretation of results.

Accessories associated with the electrodes / Compatible combinations

Only products listed in table 2 are compatible. The use of other products could be dangerous for the patient and Alcis could not be held as responsible. The verification of the compatibility of products used is under practitioner responsibility.

Cleaning, Disinfection and Sterilization

Given the nature of the act (surgery affecting the central nervous system), do respect the laws in force for the prevention of infectious risks include transmission of prions.

- **Single use devices** *any devices identified with symbol shown here*  (refer to table 2 or product label).
 - Single-use and provided sterile by ethylene oxide.
 - Do not re-use (re-use may expose patients to unacceptable risks of infection)
 - Do not re-sterilize (re-sterilization may induce damages and performance decrease)
- **Reusable surgical accessories** *any devices identified with any of the symbols shown here*   (refer to table 2 or product label).
 - Reusable devices provided unsterile
 - Before use, clean and sterilize according to a controlled and validated process (alkaline detergent and sterilization in an autoclave at 134 ° C for 18 minutes).
 - During the cleaning, take care to biological residues absence.
 - After use, clean the reusable ancillaries (alkaline detergent).
- **Extension cable range 2069-CE**
 - Reusable and provided non sterile.
 - Clean and disinfect using a cloth dampened with alcohol before and after use.

Storage and maintenance

- Store in a dry environment at room temperature.
- Inspect packaging before use to detect any damage.
- Do not use if package is opened or damaged (risk of fall of the devices which can cause a deterioration or a breakage of the device)
- **For all products provided sterile** *any devices identified with symbol shown here*  (refer to table 2 or product label).
 - ➔ Do not use after the expiry date printed on the packaging.

Patient monitoring

To avoid any risk of incidents, the patient must be placed under monitoring during the phases of recording, stimulation and coagulation.

Materials in contact with the patient

The materials in contact with brain and cerebrospinal fluid throughout implantation are:

- Platinum-iridium (alloy 90% Platinum, 10% Iridium),
- Polyurethane
- Medical Grade Silicone,
- Titanium Alloy TiAl6V4 ISO 5832-3.

Disposal

There are infectious risks related to stopper disposal, according to WHO waste classification. To minimize the risks, the most common treatments for infectious waste are: steam treatment (autoclaving), integrated steam treatment (hybrid autoclaves or advanced steam treatment), microwave treatment, dry heat treatment, chemical treatment (internal shredding or chemical disinfectants) and incineration (combustion, pyrolysis or gasification). One of this method shall be followed to reduce the infectious risk.

According to WHO, electronic devices as extension cable shall be disinfected with common disinfectant, prior be disposed by providing it to companies specialized in electronic waste disposal

There is no particular risk related to removable link disposal according to WHO waste classification. The user shall follow local legislation for waste disposal.

5. Main complications

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

The main complications associated with the implantation of depth electrodes are:

- Infections,
- Formation of hematomas and hemorrhages,
- Permanent or transitory neurological deficits.

Complications due to the intervention techniques in stereotactic conditions.

The respect for the stereotactic conditions allows to decrease the risks of complication during the implantation of electrodes. According to the state of the art, the depth electrodes being implanted since the 50s, the complications due to the interventions are well known as well as the way to manage those potential complications.

6. Warranty limit

ALCIS products are submitted to control tests at several stages of manufacturing process to ensure compliance with specifications. ALCIS can never be held responsible for any expenses or damages directly or indirectly resulting from:

- The non-respect of instructions, recommendations and warnings of the instruction of use,
- The purchase, the use, the implantation, the explantation or replacement of depth electrodes or their accessories,
- The use of accessories or devices unsuited or defective.

Use procedure

All the ancillaries (see Table 2) necessary for the implantation of the ALCIS depth electrodes in stereotactic way are proposed by ALCIS. **Contact your local reseller.** The use of other products could be dangerous for the patient and Alcis could not be held as responsible. The verification of the compatibility of products used is under practitioner responsibility.


6.1. Operating mode

Step 1. Stereotactic positioning


The method is based on a preoperative planning. It aims to define the paths (trajectories) of the electrodes to reach areas of interest while avoiding the cerebral vessels (vessel free corridors). It is performed using complementary imaging methods such as radiography, angiography, tomodensitometry and MRI.

Step 2. Implantation

S2.1 Positioning of guidance screws


The guidance screws (see table 2) must be positioned to insure the guidance of the electrodes during their implantation according to trajectories defined during stereotactic location.

Take note of the distance between the proximal edge of the guide (or grid) and the skin, (distance D_a).

Drill the bone, directly through the skin, till its internal wall (sensory feedback). Secure the depth with a suitable depth-stopper (see Table 2). Move gradually the depth-stopper to follow the drilling progress.


To avoid deviation in case of tangential trajectories, make the first cut of bone external table by small blows and with a maximum revolution speed.

Take note of the distance between the end of the drill and the depth-stopper (D_b).

Open the dura-mater with the coagulation electrode (see Table 2). Secure penetration depth, with a suitable stopper. Place the depthstopper at a distance suitable to open the dura equal to: **Db + 5mm**.


Determine **Ls** the length of guidance screw to use by adding 10mm to the bone + skin thickness ($D_b - D_a + 10mm$).

Insert the guidance screw by square end into a long screwdriver. Insert this assembly through the grid or guide in the osseous hole screw is self-threading. **Screw the screw**.

Using short guidance screw (2023-VG-C)


Determine **Ls** the length of guidance screw to use by choosing the closest longer length to the bone + skin thickness ($D_b - D_a$).

Hold the short guidance screw with the pincer 2023-PIN at the entrance hole.

Place the long screwdriver through the grid or the guide and insert it into the guidance screw (its hexagonal part entirely enter into the screw). Without the pincer, fully screw in the osseous hole (screw is self-threading).


Secure screwing depth with a suitable depth-stopper placed at a distance equal to **Db** from the distal end of guidance screw.

Reiterate the above operations for each screw.


We advise to control the screws position (matching with planned trajectories) with a suitable technique.

S2.2 Electrodes implantation

Caution: do not remove the electrode from its protective tube before having adjusted the implantation length.

Dt distance to target / **Ds** distance to screw / **Ls** length of screw

Determine the implantation length Li (implanted part) to choose the suitable electrode

$$Li = Dt - Ds$$


The implantation length have to comply with the followings conditions:

$$Li \geq MIL$$

MIL = Minimum Implantation length

$$Li \geq UEL + 10mm$$


UEL = Useful Implantation length

Carry out a « pass hole » with a stylet

To avoid the deviation risks during the electrode introduction, carry out a « pass hole » with a stylet (see Table 2). Pass hole can only be partial but it should never exceed the depth planned for electrode.

Place the depth-stopper at a maximum length equal to the implantation length **Li** planned.


Insert the stylet through the guidance screw (do not force, must pass freely the dura if well opened) and remove the stylet.

Using short guidance screw (2023-VG-C)



Place the guide in the screw

Adjust the stylet depth stopper at a maximal length of **Li + 20 mm**

Insert and remove the stylet, remove the guide.

Adjust the seal to the implantation length L_i

Insert the electrode through the guidance screw up to the seal, then screw tight the cap.

Pay attention to avoid to move the sliding seal during this step.

Check with a slight traction that the electrode is well fixed.

Using short guidance screw (2023-VG-C)

Adjust the seal location to $L_i + 2,5$ mm


Mount the cap (provided in the package) on the electrode.

Then insert the electrode, fully screw the cap maintaining the sliding seal close to the cap (if necessary use the screwdriver adapted for the cap).

Pay attention to avoid to move the sliding seal during this step.

Check with a slight traction that the electrode is well fixed.

Maintain the electrode in a groove of the cap and put the ring to keep the electrode blocked in position.

Check the electrodes position and the absence of hemorrhage with a suitable method.


Note: The use of the graduated ruler REF 2023-TG is recommended to secure the use of depth electrodes and tools by measuring and adjusting the lengths and distances

Step 3. Recording / Stimulation

The connection of electrode to the recording/stimulation system is done with a cable designed for this purpose. The use of ALCIS cables (see Table 2) allows a simple and protected connection as well as a fast identification of the contacts via:

- The electrode and cable connectors are keyed to limit the risks of misconnection.
- The blue dots of electrode and cable connectors must be faced.
- Tip contact is #1 Contact, next is #2, ...
- The cable's bananas plugs are numbered. When connected to the electrode, the plug #1 corresponds to the tip contact of the electrode, the plug #2 to the following macro-contact...

Depending on user expectation a Removable Link REF: 2023-CE-L part of Alcis accessories range could optionally be used to increase the necessary disconnection force of cables and depth electrode.

Warning: Do not pull on tubular part of electrode neither on wire of the cable in order to disconnect them.

Step 4. Analysis of recordings for the location of pathological brain areas to treat

Step 5. (optional) Thermo coagulation of pathological brain areas

Refer to **Instruction for coagulation** (ref. 2069-EPC-IFC).

Step 6. Electrodes removal

Warning: To do without delay at the end of the thermo-coagulation procedure

Check preliminary the absence of anomalies with a suitable technique.

Unscrew the cap, then pull the electrode to remove it.

Using short guidance screw (2023-VG-C)

Remove the retaining ring and unscrew the cap.

Pull the electrode to remove it

Unlock the screw with a wrench, then unscrew totally with a screwdriver.

Close the cutaneous hole.

Reiterate the same operations for each electrode.

Check the cerebral structures with a suitable technique.

Step 7. Checking the anatomical extent of the lesion by MRI

To be carried out within 2 months after the coagulation.

Step 8. Disposal

The disposal is under user's responsibility. Local legislation must be observed to avoid any infectious risks.

6.2. Electro-stimulation principle and parameters

ALCIS depth electrodes can be used to short stimulation to provoke crisis in order to locate area and ways of epilepsy propagation.

To avoid any risk associated with electrode heating, stimulation parameters should not exceed the following limit:

$$I^2 \times \Delta t \leq 0,002 \text{ J}\Omega^{-1}$$

with: I = stimulation current, expressed in amps
 Δt = effective duration of stimulation, expressed in second
 $=$ frequency (Hz) x lasted total duration of stimulation (s) x lasted of an impulse (s)

By example:

The following stimulation parameters:

$$I = 12 \text{ mA, frequency} = 50 \text{ Hz, total duration of stimulation} = 5 \text{ s, impulse duration} = 1 \text{ ms}$$

$$I^2 \times \Delta t = (0,012)^2 \times 50 \times 5 \times 0,001 = 0,000036 \leq 0,002 \text{ J}\Omega^{-1}$$

6.3. Thermo-coagulation principle

Refer to **Instruction for coagulation** (ref. 2069-EPC-IFC).






7. Definition of labelling and packaging symbols

Table 1 – Symbols definition

	Single use device/ Do not reuse		Manufacturer's Catalog Reference		Store in a dark place		Sterilized by ethylene oxide
	Use before		Batch No		Store in a dry place		Reusable device validated for x nn use cycle (here 50).
	Do not use if damaged		Period of use After Opening (here 24 months)		Steam sterilize before use: 134°C – 18min		Sterilize to re-use validated for x nn cycles (here 50)
	Manufacturer		Warning		Read Instruction For Use		Warning, refer to Instruction For Use
	Applied part type BF		Non sterile device		International phone number		MIL Minimum Implantation Length UEL Useful Exploration Length
	Medical Device		Manufacturing date				

Table 2 – Products and accessories

Range / product	Products	Characteristics / accessories compatibility	
Range 2069-EPC 	Depth electrode Technical form 2069-EPC-FT		Use exclusively with Alcis guidance screws range 2023-VG(-C). <ul style="list-style-type: none">Electric resistance contact-connector pin $\leq 140 \Omega$Refer to label to identify:<ul style="list-style-type: none">electrode type (macro; micro / macro...).exploration part specifications.
Range 2069-CE Year of CE marking: 2008	Extension cable Technical form 2069-CE-FT		Use exclusively with ALCIS product range 2069 (depth electrodes and accessories) Connection: <ul style="list-style-type: none">Specific 2069 range connectorFemale isolated banana plugs Ø1,5mm
Range 2069-CE Year of CE marking: 2008	Extension cable micro-macro Technical form 2069-CE-FT		Use exclusively with ALCIS Depth electrodes range 2069-EPC micro-macro Connection: <ul style="list-style-type: none">Specific 2069 range connectorFemale isolated banana plugs Ø1,5mm
2023-CE-L 	Removable link 		Use exclusively with ALCIS depth electrodes range 2069-EPC and Alcis extension cable range 2069-CE
Range 2023-VG Year of CE marking: 2013	Guidance screw Technical form 2023-VG-FT		Use exclusively with ALCIS depth electrodes range 2069-EPC Screw exclusively with ALCIS screwdriver range 2023-TO <ul style="list-style-type: none">Available in length 15, 20, 25, 30 and 35mm

Range 2023-VG-C  1014	Short guidance screw <i>Technical form 2023-VG-C-FT</i>		 IFU : 2023-VG-119	Use exclusively with ALCIS depth electrodes range 2069-EPC Screw exclusively with ALCIS screwdriver range 2023-TO-C • Available in length 10, 15, 20, 25 and 30mm
Range 2023-TG  1014	Graduated ruler <i>Technical form 2023-TG-FT</i>		 IFU : 2023-TG-013	Maximal measuring length : 175 mm Compatible with ALCIS depth electrodes range 2069-EPC and installation accessories range 2023 or 2036.
Range 2023-FO  1014	Drill <i>Technical form 2023-APR-FT</i>		 IFU : 2023-AS-IFU	Cutting Ø 2,1 mm Cutting length 20 mm Body OD adapted to ID 2.5mm guides. Length on request in the limit of 250 mm.
Range 2036-FO  1014	Drill (*) (without guide nor fixed stopper) <i>Technical form 2036-AS-FT</i>		 IFU : 2036-AS-IFU	Check guiding tool compatibility by referring to drill ref. engraved on body. Cutting Ø is engraved on drill body – Ø 2.1 mm fit with 2023-VG Guidance screw range. Length on request in the limit of 250 mm.
Range 2023-EC  1014	Coagulation electrode <i>Technical form 2023-EC-FT</i>			Tension max allowed 2kV Advice set up: 30W in coagulation mode
Range 2023-ST  1014	Single use stylet <i>Technical form 2023-APR-FT</i>		 IFU : 2023-ST-IFU	External diameter OD 0,8 mm Useful length refer to label With stopper (2023-ST-AB) or without (2023-ST-SB).
Range 2023-TO  1014	Screwdriver <i>Technical form 2023-APR-FT</i>		 IFU : 2023-AS-IFU	Use 2023-TO exclusively with ALCIS standard guidance screws range 2023-VG. Use 2023-TO-C exclusively with ALCIS short guidance screws range 2023-VG-C.
Range 2023-CL  1014	Wrench <i>Technical form 2023-APR-FT</i>		 IFU : 2023-AS-IFU	Use 2023-CL (useful width 1.9mm) exclusively with ALCIS standard guidance screws range 2023-VG. Use 2023-CL-C (useful width 3.5mm) exclusively with ALCIS short guidance screws range 2023-VG-C.
2036-BU-dd-RU  1014	Stopper Ø d,d mm <i>Technical form 2023-APR-FT</i>		 IFU : 2036-AS-IFU	Adapted to instruments with an OD ≈ d,d mm Use exclusively with tools with compatible OD Check tool compatibility with dd in reference. dd = nominal diameter in tenths of millimeters
2036-TRX-dd-III  1014	Testing/spoting rod <i>Technical form 2036-AS-FT</i>		 IFU : 2036-AS-IFU	Check guiding tool compatibility before use. OD=d.dmm, length = llmm Type =X, X=1 stainless steel
2036-PCX-dd-III  1014	Conical point for bone <i>Technical form 2036-AS-FT</i>		 IFU : 2036-AS-IFU	Check guiding tool compatibility before use. OD=d.dmm, length = llmm Type =X, X=1 stainless steel

(*) (Other dimensions available on request, consult us)

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